













**Laerdal**

helping save lives

# Defibrillator Training

Defibrillator	I want to use real defibrillators on CPR-D manikins	I want to use real defibrillators in training mode	I want to train with AED Trainer
<p><b>HeartStart AED</b></p>	<p>You can use any of the following rhythm generating manikins</p> <ul style="list-style-type: none"> <li>• Resusci Anne CPR-D</li> <li>• Resusci Anne Simulator</li> <li>• ALS Simulator</li> <li>• SimMan</li> </ul>	<p>Use on any manikin</p> 	<p>Use on any manikin or with Little Anne AED Training System or AED Resusci® Anne</p>
<p><b>HeartStart FR2 +</b></p>  <p>M3840ABU M3841ABU</p>	<p>945004 Manikin Electrode Adapter.</p> 	<p>M3864A FR2 Training and Administration Pack</p>  <p>07-10900 HeartStart Training Pads M3871A FR2 Infant / Child Training Pads</p> 	<p>94005033 AED Trainer 2</p> <p>02006033 AED Little Anne Training System</p> <p>07-10900 HeartStart Training Pads (No pads on/off detection)</p> <p>945090 Training Pads Laerdal Link (to use with AED Little Anne or AED RA)</p> <p>945050 Remote control for AED Trainer 2</p> <p>945055 Programming kit for AED Trainer 2</p> 
<p><b>HeartStart FRX</b></p>  <p>861305A</p>	<p>945004 Manikin Electrode Adapter.</p> 	<p>989803139271 Training Pads II (1 pair) M5089A External manikin adapters (5 on a roll)</p> 	<p>861307ABU HeartStart FRx Trainer</p> <p>989803139271 Training Pads II (1 pair) M5089A External manikin adapters (5 on a roll) 942510 FRX Trainer Soft Pack</p>
<p><b>HeartStart HSI</b></p>  <p>M5066A</p>	<p>Not applicable.</p>	<p>Need to use external or internal (not recommended) manikin adapter.</p> <p>M5073A Adult Training Pads Cartridge M5074A Infant/child training pads cartridge M5093A Replacement adult training pads M5094A Replacement infant/child training pads</p> <p>M5089A External manikin adapters (5 on a roll)</p>	<p>M5086AABU HS-I Trainer</p> <p>Need to use external or internal (not recommended) manikin adapter.</p> <p>M5073A Adult Training Pads Cartridge M5074A Infant/child training pads cartridge M5093A Replacement adult training pads M5094A Replacement infant/child training pads</p> <p>M5089A External manikin adapters (5 on a roll)</p> 